

ENDOSCOPIC TELESCOPES

INGENIOUS

1. DEFINITIONS
2. KEY COMPONENTS FOR IMAGING SYSTEM
3. ENDOSCOPE CONSTRUCTION
4. ENDOSCOPE PARAMETERS
5. ENDOSCOPE SPECIFICATIONS
6. TYPES OF TELESCOPE
7. ENDOSCOPE FEATURES
8. REPROCESSING PROCEDURE
9. CAUTIONS
10. TAKE HOME MESSAGE

- Is a minimal invasive surgery procedure (MIS)
- Is the examination and inspection of the interior of body organs, joints or cavities through an **endoscope**
- It allows physicians to peer through the body's passageways



ENT



GYNECOLOGY



LAPAROSCOPY



ORTHOPEDICS



UROLOGY


- Is an optic instrument using fiber optics and powerful **lens** systems.
- It provides lighting and visualization of the interior body.
- Magnification of the image or angulation of the reflected light




WHY DO WE USE A TELESCOPE?

- Small diameter for small incision operations
- Compatible to surgical instruments
- First image detector
- Telescope is connected to Camera Head Unit in order to transmit light





2. KEY COMPONENTS FOR IMAGING SYSTEM



MEDICAL MONITOR



CAMERA CONTROL UNIT (CCU)




CAMERA HEAD UNIT (CHU)



LIGHT SOURCE




ENDOSCOPIC TELESCOPE



- Housing from stainless steel
- Junction laserwelded and soldered

- DOV is the angle between the optical axis of the eye and optical axis of the objective/ endoscope.
- A color-coded ring situated on the light guide connector indicates the direction of view of the endoscope.

Direction of view (DOV) DOV 0° or 30°



DOV 0°



DOV 30°





Color	Direction of View
Green	0°
Red	30°
Black	45°
Yellow	70°


- FOV Describes the visible sector of the image.

FOV:

- Standard field of view
- Wide angle field of view

Field of view (FOV) 70°





The surgeon select the appropriate diameter and length of telescope depending on the type of surgery .



Locking for Instrument connection




Locking for Light cable with different adapters



INGENIOUS medical offers the below telescopes:

- A. LAPAROSCOPES
- B. CYSTOSCOPES / HYSTEROSCOPES
- C. NEPHROSCOPES
- D. URETER-RENOSCOPES
- E. ARTHROSCOPES
- F. SINUSCOPES



- Used for laparoscopic diagnostic and surgical procedures
- A special emphasis was placed on creating a homogenous and sharp-edged illumination
- Together with the optimized optical system, the result is a very bright image with rich contrast
- The system was designed so that the illumination fibers are fastened strain-free which leads to an improved lifetime.



LAPAROSCOPES diameter 2.0 mm, working length 258 mm


- Direction of view 0°, wide angle
- Direction of view 30°, wide angle

LAPAROSCOPES diameter 5.0 mm, working length 312 mm

- Direction of view 0°, wide angle
- Direction of view 30°, wide angle

LAPAROSCOPES diameter 10.0 mm, working length 344 mm

- Direction of view 0°, wide angle
- Direction of view 30°, wide angle
- Direction of view 45°, wide angle



- Used for urology and hysteroscopy procedures
- A special emphasis was placed on creating a homogenous and sharp-edged illumination
- Together with the optimized optical system, the result is a very bright image with rich contrast
- The system was designed so that the illumination fibers are fastened strain-free which leads to an improved lifetime.



CYSTOSCOPES / HYSTEROSCOPES diameter 2.7 mm, working length 205 mm



- Direction of view 0°, Standard
- Direction of view 30°, Standard

CYSTOSCOPES / HYSTEROSCOPES diameter 2.9 mm, working length 302 mm

- Direction of view 0°, wide angle
- Direction of view 30°, wide angle

CYSTOSCOPES / HYSTEROSCOPES diameter 4.0 mm, working length 302 mm

- Direction of view 0°, Standard / wide angle
- Direction of view 12°, Standard angle
- Direction of view 30°, Standard / wide angle
- Direction of view 70°, wide angle



- URS is constructed for smooth and gentle insertion
- of the special ocular construction, there is no need for camera re-focusing and no loss of picture
- Offered with detachable instrument channel
- Standard and mini URS are available





STANDARD URS

- Working diameter 7.5 / 9.8 FR, working length 310 mm
- Working diameter 7.5 / 9.8 FR, working length 425 mm

MINI URS

- Working diameter 6 / 7.5 FR, working length 310 mm
- Working diameter 6 / 7.5 FR, working length 425 mm



- Used for percutaneous nephroscopy
- Offered with detachable instrument channel
- Compatible with 18.5 Fr. and 24 Fr. nephroscopy sheaths




STANDARD NEPHROSCOPES

- Working diameter 19 FR, working length 220 mm

MINI NEPHROSCOPES

- Working diameter 13.5 FR, working length 220 mm
- Working diameter 13.5 FR, working length 310 mm



- Used for small joints or for knee, shoulder and hip
- High image quality, robustness, and state-of-the-art construction




Arthroscope diameter 2.7 mm, working length 187 mm

- Direction of view 0°, Standard
- Direction of view 30°, Standard
- Direction of view 70°, Wide angle

Arthroscope diameter 4.0 mm, working length 175 mm

- Direction of view 0°, wide angle
- Direction of view 30°, wide angle
- Direction of view 45°, wide angle
- Direction of view 70°, wide angle



- Used to examine nasal and sinus passages.
- High image quality, robustness, and state-of-the-art construction.



SINUSCOPES diameter 2.7 mm, working length 110 mm

- Direction of view 0°, Standard
- Direction of view 30°, Standard
- Direction of view 70°, Standard

SINUSCOPES diameter 4.0 mm, working length 175 mm

- Direction of view 0°, wide angle
- Direction of view 30°, wide angle
- Direction of view 30°, Standard
- Direction of view 70°, wide angle



GENERAL REPROCESSING PROCEDURE:

1. Pre-cleaning immediately after use
2. Cleaning (manual or automated)
3. Sterilization

Pre-Vacuum Process	Regular Cycles	WHO Cycle
Temperature (°C/°F)	134/273	134/273
Holding time (min)	4	18
Drying time	30	30

 AUTOMATED CLEANING AND PRE-VACUUM STEAM STERILIZATION CONFIRMS THE MOST RELIABLE RESULTS



Don't use pliers to remove adapter

Brush thread carefully



Brush thread carefully



Remove light connection adapters with rubber and clean surface with bud

Don't use sharp edge instrument for cleaning



Only use of cotton bud or tissue and alcohol for ALL cleaning

10. TAKE HOME MESSAGE - TELESCOPES

LIFE CYCLE

ILLUMINATION FIBERS ARE FASTENED STRAIN FREE. NO TENSILE STRESS WHILE AUTOCLAVATION



MINIMAL REFLECTION

BLACK MATT ROUGH SURFACE MINIMIZE THE LIGHT REFLECTION



ILLUMINATION

HOMOGENOUS AND SHARP-EDGED ILLUMINATION



OPTICS

INCREASED DEPTH OF FIELD HELPS IN CHANGING WORKING DISTANCE W/T REFOCUSING



BEST IMAGE QUALITY

SUPERIOR LENS MATERIAL



PATIENT SAFETY

SMOOTH EDGE TO AVOID TISSUE DAMAGE




UNIVERSAL CONNECTIVITY

COMPATIBLE TO ALL CAMERA HEADS FROM ANY SUPPLIER



VARIETY OF MODELS

DIFFERENT DIAMETERS & WORKING LENGTHS



MULTI - STERILIZATION MODES

HIGH & LOW TEMPERATURE STERILIZATION PROCEDURE



INGENUITY
FOR
HEALTH



INGENIOUS

THANK YOU

INGENIOUS