


INGENIOUS



ENDOSTAP

- Full Range of Reloads for:
 - Vascular / Thin Tissue
 - Medium / Thick Tissue
 - Extra Thick Tissue
- Provides a smooth firing force through tissue due to stronger knife bar
- The strength of Anvil raised 5 times after being processed with CNC technique
- Different Endo Cutter Sizes:
 - Small / 75 ± 10 mm
 - Medium / 155 ± 10 mm
 - Large / 255 ± 10 mm



INGENIOUS medical GmbH
79106 Freiburg
Germany

+49 761 15521831

www.ingenious-medical.com

info@ingenious-medical.com

INGENUITY
FOR
HEALTH